# **Heat**Sensor

HEATSENSOR® DUO - OCTO are incubators ensuring all Unisensor diagnostic tests to receive a fixed incubation temperature no matter what the ambient / operating temperature.

Both incubators are a **One Step Test Procedure**.

## **KEY-POINTS**

- → Automatic release of dipstick.
- → Compatibility with all test kits from UNISENSOR.
- → DUO is for maximum 2 tests at a time.
- → OCTO is for maximum 8 tests at a time.
- → Self-timing with light and audio signal.
- Countdown timer display.
- → Product's specific program included.
- → Easy to use.



### PRODUCTS REFERENCES

APP032

HEATSENSOR DUO

APP052

HEATSENSOR OCTO



SPECIFIC PROGRAM

Select the right program and place the dipstick.



G Heat Sensor

#### **AUTOMATIC**

The dipstick falls down automatically into the microwell after the first incubation time.



### Unisensor S.A.

Liège Science Park Allée de la Cense-Rouge, 98 B-4102 Ougrée (Liège) BELGIUM

Phone +32 4 252 66 02 Fax +32 4 252 90 55 info@unisensor.be

www.unisensor.be

