

# BioKits

## Peanut Assay Kit



### Intended Use

BioKits Peanut Assay Kit is used for the detection and quantification of peanut in foods and environmental swabs.

### Product Specifications

Lower limit of detection:	0.1 ppm peanut
Range of quantitation:	1–20 ppm
Testing time:	75 minutes
Wells per kit:	96

### Results



### Ordering Information

Prod#.	Product description
902048Q	BioKits Peanut Assay Kit
901042J	BioKits Allergen Swabbing Kit – Contains 20 swabs, tubes and buffer

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)