

**Reveal® 2.0**

For *E. coli* O157:H7



**Intended Use**

The Reveal 2.0 for *E. coli* O157:H7 test provides for the rapid recovery of *Escherichia coli* O157:H7 and *Escherichia coli* O157:NM organisms in beef matrices facilitating detection and presumptive identification of the target organisms in as little as 12 hours.

**Product Specifications**

- Sensitivity: 1 cfu/analytical unit  
10<sup>4</sup> cfu/mL post enrichment
- Testing time: 15 minutes
- Tests per kit: 20
- Kit contains: 20 Reveal 2.0 for *E. coli* O157:H7 test devices  
20 disposable transfer pipettes  
20 sample cups  
1 promoter reagent

**Results**

Positive



Negative

**Ordering Information**

Prod.#	Product description
9734	Reveal 2.0 for <i>E. coli</i> O157:H7 – A one step assay providing visual results in 15 minutes after enrichment. Kit contains 20 devices.
9727	Reveal 2.0 <i>E. coli</i> O157:H7 Unitized Medium – Contains 20 bottles to enrich 20 - 65 g samples
9728	Reveal 2.0 <i>E. coli</i> O157:H7 Medium – 500 g
9475	Rack for Reveal sample cups
9738	Sample bag – 500 mL (250)